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
1. Jamison J, Maguire S, McCann J. Catheter policies for management of long term voiding problems in adults with neurogenic bladder disorders. *Cochrane Database Syst Rev* 2013;(11):CD004375.
2. Pannek J, Blok B, CastroDiaz D, et al. EAU Guidelines on neurogenic lower urinary tract dysfunction, 2013. Available at: [http://uroweb.org/wpcontent/uploads/20\\_NeurogenicLUTD\\_LR.pdf](http://uroweb.org/wpcontent/uploads/20_NeurogenicLUTD_LR.pdf) [Accessed 13 April 2016].
3. DeVivo MJ, Krause JS, Lammertse DP. Recent trends in mortality and causes of death among persons with spinal cord injury. *Arch Phys Med Rehabil* 1999; 80(11):14119.
4. Whiteneck GG, Charlifue SW, Frankel HL, et al. Mortality, morbidity, and psychosocial outcomes of persons spinal cord injured more than 20 years ago. *Paraplegia* 1992;30(9):61730.
5. Soden RJ, Walsh J, Middleton JW, et al. Causes of death after spinal cord injury. *Spinal Cord* 2000; 38(10):60410.
6. Cameron AP, Wallner LP, Tate DG, et al. Bladder management after spinal cord injury in the United States 1972 to 2005. *J Urol* 2010;184(1):2137.
7. Weld KJ, Dmochowski RR. Effect of bladder management on urological complications in spinal cord injured patients. *J Urol* 2000;163(3):76872.
8. Shavelle RM, Paculdo DR, Tran LM, et al. Mobility, continence, and life expectancy in persons with Asia Impairment Scale Grade D spinal cord injuries. *Am J Phys Med Rehabil* 2015;94(3):18091.
9. Wyndaele JJ, Maes D. Clean intermittent selfcatheterization: a 12year follow up. *J Urol* 1990; 143(5):9068.
10. PerrouinVerbe B, Labat JJ, Richard I, et al. Clean intermittent catheterisation from the acute period in spinal cord injury patients. Long term evaluation of urethral and genital tolerance *Paraplegia* 1995;33(11):61924.
11. Wyndaele JJ. Complications of intermittent catheterization: their prevention and treatment. *Spinal Cord* 2002; 40(10):53641.
12. BieringSørensen F, Nielans HM, Dørflinger T, et al. Urological situation five years after spinal cord injury. *Scand J Urol Nephrol* 1999;33(3):15761.
13. Krebs J, Wöllner J, Pannek J. Risk factors for symptomatic urinary tract infections in individuals with chronic neurogenic lower urinary tract dysfunction. *Spinal Cord* 2015 Dec 1. doi: 10.1038/sc.2015.214. [Epub ahead of print].
14. Esclarín De Ruz A, García Leoni E, Herruzo Cabrera R. Epidemiology and risk factors for urinary tract infection in patients with spinal cord injury. *J Urol* 2000;164(4):128589.
15. Wyndaele JJ, De Sy WA, Claessens H. Evaluation of different methods of bladder drainage used in the early care of spinal cord injury patients. *Paraplegia* 1985;23(1):1826.
16. McGuire EJ, Savastano J. Comparative urological outcome in women with spinal cord injury. *J Urol* 1986; 135(4):7301.
17. Turi MH, Hanif S, Fasih Q, et al. Proportion of complications in patients practicing clean intermittent selfcatheterization (CISC) vs indwelling catheter. *J Pak Med Assoc* 2006;56(9):4014.
18. Shekelle PG, Morton SC, Clark KA, et al. Systematic review of risk factors for urinary tract infection in adults with spinal cord dysfunction. *J Spinal Cord Med* 1999;22(4):25872.
19. Giannantoni A, Scivoletto G, Di Stasi SM, et al. Clean intermittent catheterization and prevention of renal disease in spinal cord injury patients. *Spinal Cord* 1998;36(1):2932.
20. Killorin W, Gray M, Bennett JK, et al. The value of urodynamics and bladder management in predicting upper urinary tract complications in male spinal cord injury patients. *Paraplegia* 1992; 30 (6):43741.
21. Bakke A, Vollset SE. Risk factors for bacteriuria and clinical urinary tract infection in patients treated with clean intermittent catheterization. *J Urol* 1993;149(3):52731.
22. Bakke A, Digranes A, Hoisaeter PA. Physical predictors of infection in patients treated with clean intermittent catheterization: a prospective 7 year study. *Br J Urol* 1997;79(1):8590.
23. Vahr S, CobussenBoekhorst H, Eikenboom J, et al. Evidencebased guidelines for best practice in urological health care. Catheterisation; urethral intermittent in adults; dilatation, urethral intermittent in adults. European Association of Urology Nurses (EAUN) 2013. Available at: [http://patients.uroweb.org/wpcontent/uploads/CatheterisationUrethralIntermittentinadults\\_Lr\\_DEF.pdf](http://patients.uroweb.org/wpcontent/uploads/CatheterisationUrethralIntermittentinadults_Lr_DEF.pdf) [Accessed 13 April 2016].
24. Zowawi HM, Harris PN, Roberts MJ, et al. The emerging threat of multidrugresistant Gramnegative bacteria in urology. *Nat Rev Urol* 2015;12(10):57084.
25. World Health Organization. Antimicrobial resistance Fact sheet N° 194, updated April 2015. Available at: <http://www.who.int/mediacentre/factsheets/fs194/en/> [Accessed 13 April 2016].
26. Consortium for Spinal Cord Medicine. Bladder management for adults with spinal cord injury: a clinical practice guideline for healthcare providers. *J Spinal Cord Med* 2006; 29(5):52773. Available at: <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC1949036/pdf/i10790268295527.pdf> [Accessed 25 April

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The Coloplast logo is a registered trademark of Coloplast A/S. © 2017-6 All rights reserved Coloplast A/S